

TEST REPORT

(Electronic version)

扫码下载报告



VERIFICATION WEBSITE: www.gtgc.net.cn

VERIFICATION CODE: JUVF-8061-54

No: 210089916

ISSUE DATE: 2021-05-11

APPLICANT: WUJIANG TUTAIKE TEXTILES & FINISHING CO., LTD
ADDRESS: NO. 1599, 3RD SOUTHERN RING ROAD, SHENGZE TOWN, WUJIANG, JIANGSU, 215200, CHINA

INFORMATION CONFIRMED BY APPLICANT:

ISOLATION GOWN

QUANTITY: 1 PIECE

STANDARD ADOPTED:

YY/T 0506.2-2016 <Surgical drapes, gowns and clean air suits for patients, clinical staff and equipment-Part 2:Performance requirements and test methods>

DATE RECEIVED/DATE TEST STARTED: 2021-04-27

CONCLUSION:

Cleanliness-microorganism M

NOTE: "M" -MEET THE STANDARD'S REQUIREMENT "F" -FAIL TO MEET THE STANDARD'S REQUIREMENT
"---" -NO COMMENT

REMARK:

ALL THE TESTED ITEMS ARE TESTED UNDER THE STANDARD CONDITION (EXCEPT FOR INDICATION).
COPIES OF THE REPORT ARE VALID ONLY RE-STAMPED.
THE EXPERIMENT WAS CARRIED OUT AT No. 1, ZHUJIANG ROAD, PANYU DISTRICT, GUANGZHOU, GUANGDONG, P. R. CHINA.

APPROVED BY:

Nan Ma ENGINEER

马楠



TEST REPORT

(Electronic version)

No:210089916



TEST REPORT

(Electronic version)

No:210089916

Cleanliness-microorganism

M

(YY/T 0506.7-2014, Eluting parameter: Frequency: 200r/min, Time: 1min)
Requirement

Total number of bacterial colony eluted:
0 CFU

Recovery rate:
80%

Correction factor:
1.25

Cleanliness-microorganism:
0 CFU/dm²

≤300 CFU/dm²
(Standard performance)
(YY/T 0506.2-2016)

—End of Report—

PAGE 3 OF 3

